






Capturing and Analyzing Bragg Peak Proton Range Data on the Mevion HyperScan

XRV-3000 Eagle



using the XRV-3000 Eagle and LCW-300 multi-chevron wedge

XRV-3000 Eagle Setup




Tallest LCW-300 chevron should be located at the rear of the Eagle facing the handle

The LCW is centered on the target and the Eagle is aligned so the lasers hit the chevron peaks.

Mevion LCW-300 Delivery Plan

> 190 MeV,
9x12 - 15mm pitch



The grid pattern for each energy layer approximates a flat field


Dose for each spot is ~0.2 MU

There should be no beam activity in the capture ROI for ~1 second so the Logos software can detect the change in beam energy


The standard spreadsheet analysis template allows for 35 energy levels

Mevion LCW-300 Delivery Plan


< 165 MeV,
6x8 - 25mm pitch




165 – 190 MeV,
7x9 - 20mm pitch




> 190 MeV,
9x12 - 15mm pitch




Each energy layer is captured as a proton radiograph allowing 2 peak-to-peak measurements





Pristine Bragg Peak

Quenched Bragg Peak




The pristine Bragg Peaks are “quenched” in the scintillator causing the peak-to-peak measurements to be skewed in the proximal direction.

It has been experimentally observed that acrylic and Teflon LCW measurements when adjusted for the WET, correspond closely to the proximal 80% range of protons in water.

Setting up BeamWorks Strata

Launch BeamWorks Strata

Select “Custom Scripts”




Capture Script Style:

- Generate Beamcenter (Automatic)
- Capture Images to Working Data
- Process Images from Working Data
- Real-Time Memory Capture
- Capture Peak Image [Manual]
- Custom Scripts

Select “Browse”


Setting up BeamWorks Strata

Load BraggPeakCaptureEagle
script



Select "Open"



Setting up BeamWorks Strata



Set the Eagle camera
Frame Rate to 2 frames
per second and the Gain
to 24 dB



Making Measurements

Click on Start Script to launch WinLVS




Making Measurements

WinLVS captures images and saves measurements to disk




Tap spacebar to add information windows


Deliver plan



Click "Option 1"



Making Measurements



Pressing the Option 1 button starts the process of detecting Bragg Peak regions in the captured images and writing peak-to-peak measurements to file in a CSV format that can be easily imported into the template spreadsheet for analysis.

Data is then pasted into the analysis spreadsheet in order to calculate the equivalent proton range in water

	A	B	C	D	E	F	G	H	I	J	K
10		Folder 0008		Layer Energy (MeV): 170.9 Snout Adj							
11	1	LCW-200/300R		Center	Peak-to-Peak		Chevron Height	Target	Scint. Cov	Wedge	
12	Image:	4		Dist (mm)	Dist (mm)		(mm)	(mm)	(mm)	WET	
13	Chevron: A			0	0		170.2		3.4	1.8	1.821
14	Chevron: B			0	0		150		3.4	1.8	1.821
15	Chevron: C			-12.462	62.3974		130		3.4	1.8	1.821
16	Chevron: D			14.1868	16.9677		110.1		3.4	1.8	1.821
17	Chevron: E			0	0		90.1		3.4	1.8	1.821
18	Chevron: F			0	0		70.2		3.4	1.8	1.821
19											
20		Folder 0006		Delta (MeV) 0 Layer Energy (MeV): 170.9 Snout Adj							
21	2	LCW-200/300I		Center	Peak-to-Peak		Chevron Height	Target	Scint. Cov	Wedge	
22	Image:	1		Dist (mm)	Dist (mm)		(mm)	(mm)	(mm)	WET	
23	Chevron: A			0	0		170.2		3.4	1.8	1.821
24	Chevron: B			0	0		150		3.4	1.8	1.821
25	Chevron: C			-11.9351	60.4817		130		3.4	1.8	1.821
26	Chevron: D			13.7981	16.9677		110.1		3.4	1.8	1.821
27	Chevron: E			0	0		90.1		3.4	1.8	1.821
28	Chevron: F			0	0		70.2		3.4	1.8	1.821
29											
30		Folder 0008		Delta (MeV) 10.5 Layer Energy (MeV): 160.4 Snout Adj							
31	3	LCW-200/300I		Center	Peak-to-Peak		Chevron Height	Target	Scint. Cov	Wedge	
32	Image:	3		Dist (mm)	Dist (mm)		(mm)	(mm)	(mm)	WET	
33	Chevron: A			0	0		170.2		3.4	1.8	1.821
34	Chevron: B			0	0		150		3.4	1.8	1.821
35	Chevron: C			-11.2606	82.3755		130		3.4	1.8	1.821
36	Chevron: D			14.9266	40.5036		110.1		3.4	1.8	1.821
37	Chevron: E			0	0		90.1		3.4	1.8	1.821
38	Chevron: F			0	0		70.2		3.4	1.8	1.821

Parameters for WET (Water Equivalent Thickness)




entered into the spreadsheet:

1.821 for Teflon (LCW-300 Chevrons)

1.172 for Acrylic (Eagle support Target)

K	L	M	N	O	P	Q	R	S
		Global Input Values						
		XRV-3000 SAD X mm:		1850				
		XRV-3000 SAD Y mm:		1850				
		Snout Adj. MeV:		3.1				
		Nominal Wedge WET:		1.821				
		LCW-300 Chevrons C and D.		Nominal TARGET WET				
				1.172				
		approximately 20 mm more air to reach each lower Chevron.						
		Snout Adj	3.1	Final MeV	167.8	Bragg Depth (mm):	PSI Equation for Proximal 80% F	
		Wedge	PMMA	Trap. Scaling	Bragg Penetration		Adjusted BPD	
		WET	WET	with SAD X	Depth in mm		(mm)	
		1.821	1.172	0.939081				
		1.821	1.172	0.95				
		1.821	1.172	0.960811		188.238		188.209
		1.821	1.172	0.971568		191.5767		191.5691
		1.821	1.172	0.982378				
		1.821	1.172	0.993135				

Optimization of WET values with the “Hill Climbing” method




Bad trial

Good trial

Correct WET

Excellent agreement between Mevion and Varian energy layers



Delta (MeV)	2.3	Layer Energy (MeV):	139.1	Snout Adj	2.1	Final MeV	137 Br	Bragg Depth (mm):	PSI Equation for Proximal 80% Penetration Depth:	133.86
Peak	Chevron Height (mm)	Target (mm)	Scint. Cov (mm)	Wedge WET	PMMA WET	Trap. Scaling with SADY		Bragg Penetration Depth in mm		Adjusted BPD (mm)
	170.7	3.4	1.8	1.821	1.172	0.946095				
	150.5	3.4	1.8	1.821	1.172	0.955714				
	130.4	3.4	1.8	1.821	1.172	0.965286				
	110.3	3.4	1.8	1.821	1.172	0.974857		133.8305		133.8284
	90.2	3.4	1.8	1.821	1.172	0.984429		134.2768		134.2529
	70.1	3.4	1.8	1.821	1.172	0.994				
							Average	134.0537	Average	134.0406

Delta (MeV)	10.4	Layer Energy (MeV):	140.1	Snout Adj	3.1	Final MeV	137 Br	Bragg Depth (mm):	PSI Equation for Proximal 80% Penetration Depth:	133.86
Peak	Chevron Height (mm)	Target (mm)	Scint. Cov (mm)	Wedge WET	PMMA WET	Trap. Scaling with SAD X		Bragg Penetration Depth in mm		Adjusted BPD (mm)
	170.2	3.4	1.8	1.821	1.172	0.939081				
	150	3.4	1.8	1.821	1.172	0.95				
	130	3.4	1.8	1.821	1.172	0.968811				
	110.1	3.4	1.8	1.821	1.172	0.978568		133.9582		133.9232
	90.1	3.4	1.8	1.821	1.172	0.982378		134.1826		134.1458
	70.2	3.4	1.8	1.821	1.172	0.993135				
							Average	134.0704	Average	134.0345

Two different LCW units,
 Two different XRV phantoms,
 Two different proton delivery systems,
 Two different facilities.
 Excellent agreement!

Calculation for depth of 80% energy loss of beam in water converts MeV (black) to millimeters (red)

M	N	O	P	Q	R	S	T	U	V	W	
Global Input Values											
XRV-3000 SAD X mm:		1850								The plan layer energies were de	
XRV-3000 SAD Y mm:		1850								The delivery was a grid pattern	
Snout Adj. MeV:		3.1								The PSI R80 depth equation cor	
Nominal Wedge WET:		1.821								Proximal edge R80 is used beca	
Nominal TARGET WET		1.172									
re air to reach each lower Chevron.											
Final MeV		167.8	Bragg Depth (mm): PSI Equation for Proximal 80% Penetration Depth:								190.8867
Trap. Scaling			Bragg Penetration			Adjusted BPD					
with SAD X			Depth in mm			(mm)				Y (mm)	
0.939081										0	
0.95										0	
0.960811			188.238			188.209				29.97605	
0.971568			191.5767			191.5691				8.242634	
0.982378										0	
0.993135										0	
	Average		189.9073		Average	189.8891					
Final MeV		167.8	Bragg Depth (mm): PSI Equation for Proximal 80% Penetration Depth:								190.8867
Trap. Scaling			Bragg Penetration			Adjusted BPD					
with SAD X			Depth in mm			(mm)				Y (mm)	
0.939081										0	
0.95										0	
0.960811			189.9139			189.8865				29.05574	
0.971568			191.5767			191.5694				8.242634	
0.982378										0	
0.993135										0	

PSI Equation

167.8

190.8867

Calculate proton range in water from the LCW measurements using WET and other corrections

WET*[(Scaling ratio)*(Peak to peak)/2)+(Target thickness)]

M	N	O	P	Q	R	S	T	U	V	W	
Global Input Values											
XRV-3000 SAD X mm:		1850								The plan layer energies were del	
XRV-3000 SAD Y mm:		1850								The delivery was a grid pattern o	
Snout Adj. MeV:		3.1								The PSI R80 depth equation con	
Nominal Wedge WET:		1.821								Proximal edge R80 is used becau	
Nominal TARGET WET		1.172									
to reach each lower Chevron.											
Final MeV	167.8	Bragg Depth (mm):	PSI Equation for Proximal 80% Penetration Depth:					190.887			
Trap. Scaling with SAD X		Bragg Penetration Depth in mm		Adjusted BPD (mm)						Y (mm)	
0.93908										0	
0.95										0	
0.96081		188.238		188.209						29.976	
0.97157		191.577		191.569						8.24263	
0.98238										0	
0.99314										0	
Average		189.907		Average	189.889						
Final MeV	167.8	Bragg Depth (mm):	PSI Equation for Proximal 80% Penetration Depth:					190.887			
Trap. Scaling with SAD X		Bragg Penetration Depth in mm		Adjusted BPD (mm)						Y (mm)	
0.93908										0	
0.95										0	
0.96081		189.914		189.887						29.0557	
0.97157		191.577		191.569						8.24263	
0.98238										0	
0.99314										0	
Average		190.745		Average	190.728						

An automated treatment plan allows fast analysis of 35 energy layers as a proton range constancy check with an accuracy of about 0.5 mm


The screenshot displays an Excel spreadsheet titled "LCW300ProtonBeamBraggPeakAnalysis_HPTC_FixedSnoutMeV/Drop_021_2018 - Excl". The spreadsheet is organized into sections for different proton beam energies, with a summary table on the right side circled in red.

Summary Table: LCW-300 Chevron Data Summary

Plan Layer	Post Snout Energy (MeV)	Post Snout Energy (MeV) Data	Chevron BFD (mm) Data	Avg Energy BFD (mm) Data	Delta	Absolute Delta			
1	889	106.05	228.95	230.411	1.563	1.563			
2	871.1	104.95	-13*	225.3	-3.95*	228.308	-4.085	10281	10281
3	852.2	103.05	-13*	221.4	-3.001*	222.273	-4.053	0.8759	0.8759
4	833.3	101.15	-13*	217.26	-4.136*	218.251	-4.022	0.9504	0.9504
5	814	100.25	-13*	213.46	-3.701*	214.261	-3.59	0.7611	0.7611
6	795.5	100.35	-13*	209.6	-3.882*	210.302	-3.959	0.7044	0.7044
7	777.5	100.35	-2*	205.56	-4.006*	206.17	-4.33	0.5876	0.5876
8	759.5	100.35	-19*	201.57	-4.017*	202.276	-3.894	0.711	0.711
9	741.5	100.35	-2*	197.67	-3.897*	198.27	-4.064	0.5443	0.5443
10	723.5	100.35	-2*	193.32	-3.348*	194.36	-4.029	-0.136	0.1362
11	705.5	100.35	-2*	189.34	-3.385*	190.391	-3.553	-0.748	0.748
12	687.5	100.35	-2*	185.34	-4.253*	186.233	-3.957	0.1475	0.1475
13	669.5	100.35	-2*	181.35	-3.338*	182.371	-3.522	0.764	0.764
14	651.5	100.35	-2.1*	177.21	-3.333*	178.232	-4.079	0.0233	0.0233
15	633.5	100.35	-2*	173.27	-3.338*	174.386	-3.848	0.7528	0.7528
16	615.5	100.35	-16*	169.23	-2.903*	171.337	-3.053	0.043	0.043
17	597.5	100.35	-16*	165.07	-3.223*	168.302	-3.029	0.2363	0.2363
18	579.5	100.35	-16*	161.04	-2.626*	165.295	-3.006	-0.744	0.7444
19	561.5	100.35	-16*	157.08	-3.461*	162.312	-2.983	0.3339	0.3339
20	543.5	100.35	-16*	153.25	-2.624*	159.363	-2.96	-0.001	0.0076
21	525.5	100.35	-16*	149.46	-2.934*	156.436	-2.536	-0.744	0.7437
22	507.5	100.35	-16*	145.73	-3.132*	153.523	-2.919	-0.524	0.5241
23	489.5	100.35	-2.2*	142.05	-3.977*	150.536	-3.967	-0.534	0.5345
24	471.5	100.35	-13*	138.36	-3.495*	147.546	-3.39	-0.41	0.41
25	453.5	100.35	-2.6*	134.65	-4.639*	144.561	-4.555	-0.296	0.2963
26	435.5	100.35	-2.2*	131.08	-3.895*	141.571	-3.991	-0.271	0.2716
27	417.5	100.35	-2.3*	127.53	-3.977*	138.575	-3.956	-0.251	0.2509
28	399.5	100.35	-2.3*	124.01	-3.696*	135.588	-3.507	-0.302	0.3022
29	381.5	100.35	-2.4*	120.51	-3.895*	132.544	-4.024	-0.47	0.4701
30	363.5	100.35	-2.3*	117.03	-3.495*	129.039	-3.805	-0.82	0.8204
31	345.5	100.35	-2.4*	113.57	-3.895*	125.522	-3.917	-0.386	0.3864
32	327.5	100.35	-2.5*	110.14	-4.04*	122.041	-4.022	-0.594	0.5941
33	309.5	100.35	-2.4*	106.73	-4.125*	118.596	-3.956	-0.584	0.5837
34	291.5	100.35	-2.5*	103.34	-3.632*	115.187	-3.904	-0.576	0.5761
35	273.5	100.35	-1.3*	100.0	-2.088*	111.815	-2.006	-0.534	0.5342
Average									0.4947
Std Dev.									0.3534
Min.									1.003
Max.									0.0095

Total Depth: 94.822 mm Average Layer Depth: 3.647 mm

Digital Real-Time X-ray and Proton Beam Metrology Solutions



Logos
Systems

www.logosvisionsystem.com